



Atopy South East Asia Elite

Efficient allergy diagnostics with a customised multiparameter test

Blot	Code	Allergen name
Atopy South East Asia Elite	f1	Egg white
	f75	Egg yolk
	f3	Codfish
	f23	Crab
	f2	Cow's milk
	f76	nBos d 4 alpha-Lactalbumin (Milk)
	f77	nBos d 5 beta-Lactoglobulin (Milk)
	f78	nBos d 8 Casein (Milk)
	e204	nBos d 6 BSA (Milk)
	f4	Wheat flour
	f5	Rye flour
	f10	Sesame
	f13	Peanut
	f14	Soybean
	f17	Hazelnut
	f20	Almond
	f256	Walnut
	f25	Tomato
	f31	Carrot
f35	Potato	
f85	Celery	
f33	Orange	
f49	Apple	
e1	Cat	
e2	Dog	
e3	Horse	
e6	Guinea pig	
e82	Rabbit	
e84	Hamster	
m1	Penicillium notatum	
m2	Cladosporium herbarum	
m3	Aspergillus fumigatus	
m6	Alternaria alternata	
d1	Dermatophagoides pteronyssinus	
d2	Dermatophagoides farinae	
g1	Sweet vernal grass	
g2	Bermuda grass	
g3	Orchard grass	
g4	Meadow fescue	
g5	Rye grass	
g6	Timothy grass	
g12	Cultivated rye	
g13	Velvet grass	
g14	Oat	
g16	Meadow foxtail	
w6	Mugwort	
w9	Plantain (English)	
t2	Alder	
t3	Birch	
t4	Hazel	
t7	Oak	
CCD	CCD marker	
Ind	Indicator	

- **Simultaneous analysis** of specific IgE against 51 allergens
- **Small amount of serum** (100–400µl), ideal for paediatrics
- **Standardised evaluation** in EAST classes using EUROLineScan software
- **Fast & simple** test performance (< 3h)
- **Unique feature:** integrated CCD marker on each test for improved interpretation of cross reactions
- **Fully automated** incubation possible with the EUROBlotOne

